

# Progression Length Recommendation

Presio i Digital offers three different progression lengths: 10 mm\* 12 mm and 14 mm. The following guidelines may help you make the right choice for your patients:


		10 mm*	12 mm	14 mm
Patient use prior to Presio i Digital	Progressive lens ≥ to 14mm			•
	Progressive lens ≤ to 13mm		•	
	Bifocal lenses	•	•	
	Single vision distance lenses			
		10 mm*	12 mm	14 mm
Patients primary lens use	Distance vision			•
	Intermediate vision			•
	Near vision	•	•	
	Small framo			

Former reader/task-specific user

# Minimum Recommended Fitting Height

Presio i Digital is available in three progression lengths:

Progression Length	Fitting Height				
Presio i Digital 10*	13mm	$\neg$			
Presio i Digital 12	15mm				
Presio i Digital 14	17mm				



\* 10mm Progression Length Available First Quarter 2011





# Introducing Nikon's New Digital Progressive Lens: Presio i Digital with Enhanced Digital Technology

Nikon introduces Presio i Digital, a new generation of progressive lenses that integrates Nikon Advanced Digital Surfacing Technology and Optimization, to achieve a state-of-the-art lens design with increased image resolution and precision. Presio i Digital is engineered to provide wearers with sharper and wider vision.

## **Product Concept: Advanced Digital Surfacing Technology and Optimization**

Presio i Digital, Nikon's core lens within the digital line-up, follows the technological footsteps of the previous design philosophy while now offering greater precision and sharper vision. Designed to provide accurate power in the rotational wearing state, Presio i Digital utilizes optical compensation calculations to ensure the precise and correct power is processed.

## **Advanced Digital Surfacing Technology**

+ Nikon Optimization

Nikon lens designs are a result of mastering aspheric lens technology, and are based on the relationship between human visual requirements and enhanced visual comfort.



Nikon Advanced Digital Surfacing Technology encompasses the latest advancements in digital processing. Designed only once the lens order is received, the optimization process reduces as many aberrations as possible, resulting in an optimized surface and better controlled visual





## **Sharper and wider vision:**





**Conventional Progressive Lens** 






Presio i Digital

#### **Dual Indication:**

"Dual Indication" is the term used to identify that the written prescription power and the actual power read on the lens through the focimeter will be different.

- Distance vision area is increased by up to 12%
- Near vision area is increased by up to 50%
- Aberration level is reduced by up to 16%




As the prescription strengthens, the difference between prescription power in the wearing state and that measured by the focimeter increases.

This variation in power is a result of the difference in incident angle and vertex distance at the distance vision power verification point, from when the lens is in the focimeter to when the lens is in the wearing position.

Presio i Digital considers optical performance in the wearing state as a critical component of the design criteria, resulting in minimized aberrational astigmatism. Therefore we provide prescription power in the wearing state and measured power with the focimeter: dual indication.

#### **Better Control in any Prescription**

Nikon's advanced digital surfacing technology is the backbone of Presio i Digital, providing increased image resolution with more precision. This combined with Nikon's sophisticated calculation engine and the latest processing technology provides enhanced visual performance and comfort.



Digitally surfaced lenses provide more precision and control of cylinder power and axis than conventionally surfaced

lenses. This is due to the limitation that toric curvatures impose on traditional surfacing.

#### **Features & Benefits**

